

**Single Beam Microprocessor UV-VIS
Spectrophotometer LI-285
with Software**



Single Beam Microprocessor UV-VIS
Spectrophotometer LI-285
with Software

Technical Parameters:

Model	LI-285
Optical system:	CT grating monochromator
Wavelength range:	195 nm - 1020nm
Wavelength accuracy:	+/- 2 nm
Spectral bandwidth	2nm
Wavelength repeatability:	1 nm
Stray light:	<= 0.3%T @220nm
Photometric accuracy:	+/- 0.5%T
Photometric repeatability:	0.2%T
Stability:	+/- 0.004 A/h @500nm
Work Pattern:	T, A, C
Display area:	0- 200%T, -0.3 - 3A
Display modes:	4 LCD
Zero setting:	automatic
Weight:	12 kg


APPLICATIONS

- Medicine/Pharmaceutical
- Industry/Environmental Monitoring
- Commodity Inspection
- Food Inspection
- Agricultural Chemistry
- Teaching in Colleges & Universities
- Metallurgy
- Geology
- Machine Manufacturing
- Petrochemical Industries
- Water and Waste Water Labs
- Food and Beverage Labs

STANDARD CONFIGURATION

Main Set	: 1 No.
Software RS 232 Cable	: 1 No.
Accessory kit, Including a 4 glass cell	: 1 No.
2 Quartz Cell	: 1 No.
Operating Manual Software Manual	: 1 No.
Software CD	: 1 No.
Power Cable	: 1 No.

Lasany International



Mob: +91-9357947917, +91-8725830111
 Email: sales@lasanyspectrophotometers.com
 Website: www.lasanyinternational.com
 www.spectrophotometers.in
 www.spectrophotometersindia.com
 Skype ID: lasany.international