

Lasany[®]

DIGITAL VISCOMETER MODEL : LI-V201



Applications

- Dairy Products
- Juices, etc
- Pharmaceuticals
- Coatings Solvents
- Polymer Solutions Oils Paints
- Inks Latex Adhesives(Solvent base)

Lasany International



Mob: +91-9357947917, +91-8725830111

Email: sales@lasanyspectrophotometers.com

Website: www.lasanyspectrophotometers.com

www.spectrophotometers.in

www.spectrophotometersindia.com

Skype ID: lasany.international

Perfection in
Laboratory Science

CE ISO 9001 : 2015
Certified Company

Features

Stepping Motor means Accurate, reliable operation
 Direct readout of all measurement parameters
 Auto Range Showing
 Time Function for measurement
 Sound alarm at under 20% Torque
 Linear calibration
 Wide range power supply: 100V-240V



Model	LI-V201
Measuring Range(mPa.s)	1-2000K
RPM	0.3, 0.6, 1.5, 3, 6, 12, 30, 60
Spindles	#1,#2,#3,#4 is standard configuration ,(#0 is for option)
Accuracy	±1.0% of range
Repeatability	±0.5%
Display	LCD
Temperature Monitoring	Yes

Optional Accessories

- Circulating Water bath
- Low Viscosity Adapter
- Standard Oils
- Micro Printer
- RTD Temperature Probe
- Glass Beaker With Circulation

K=1 Thousand